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# Multisil-Malfarben

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## Instructions for use

# Multisil stains

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

### 1. Description of the product

Multisil stains based on wax are available in 10 different shades: white, yellow, red, blue, black, green, purple, light brown, brown and dark brown.

### 2. Indication range

Multisil stains are used for colour surface characterization of epitheses.

### 3. Contraindication

Based on the current state of knowledge there are no known contraindications.

### 4. Safety recommendations and hazard warning

After skin contact wash with copious amount of water and soap. After eye contact, rinse with copious amount of water. Please request the material safety data sheets separately.

### 5. Storage and durability

#### 5.1 Storage

Multisil stains must be stored under dry and cool conditions; the container (bottle) must be sealed after the use.

#### 5.2 Durability

In the original, sealed condition, see label.

### 6. Processing

Apply Multisil stains onto the degreased surface of the epithesis and spread it until the desired intensity is obtained. Once characterization of the epithesis is completed, apply Multisil sealing agent (REF 520 0100 5) or matt sealing agent for epithetics (REF 540 0109 4) onto the Multisil stain to prevent the stain from coming off.

### 7. Symbols

Can be used until YYYY-MM-DD

LOT Lot number

### 8. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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